Below is a list of terms relevant to the

∞ MetaMoreFore ∞

<u>Skin Tube</u>: This is the column, or tube, of skin that is created either by pulling the slack shaft skin forward of the glans or tensioned with a foreskin restoration device. It is similar to the foreskin in the intact state.

<u>Cylinder of rings</u>: When installed, the multiple rings of the *MetaMoreFore* create a cylinder within which the skin tube resides.

<u>Bottom and Top End Rings</u>: These are components and form the "ends" of the cylinder of rings as created with the **MetaMoreFore**.

<u>Set Rings</u>: These are components of the **MetaMoreFore** and are slip over the skin tube as a means of constructing the cylinder of rings. *Set Rings* come in five styles: *Standard*, *Wedged*, *DILE*, *Weighted*, and *Extreme Tensioning*.

<u>Tensioning Plug</u>: This component of the **MetaMoreFore** is a slightly conical shaped plug that serves to keep skin tensioned. In other restoration devices, the "gripper" serves this purpose. There is a hole in the *Tensioning Plug* to allow for the evacuation of urine.

<u>Set</u>: This refers to the construction of the initial cylinder of rings, which includes the *Bottom End Ring*, any installed *Set Rings* and the *Top End Ring*.

<u>Set Length</u>: This refers to the full length of all *Set Rings* that have been installed during the "set." This length does not include the length of the *Bottom* and *Top End Rings* (which in composite account for ~1.8 cm of the full length upon "set"), but for convenience and since all "sets" require these rings, their length is not added to the "set length." For instance, a "set" using a 2.0 cm as well as a 1.0 cm *Set Ring* yields a "set length" of 3.0 cm. But in actuality, ~5.8 cm of a skin tube has been created (~1.8 cm from the *End Rings* and 3.0 cm from the *Set Rings*).

<u>Expansion Rings</u>: These are components of the **MetaMoreFore** and are inserted onto the skin tube (between two *Set Rings*) once the **MetaMoreFore** has been "set." These rings help to further tension the **MetaMoreFore**. Expansion Rings come in three styles: Standard, Wedged, Weighted and Strapping-Specific (Standard or Wedged).

<u>Strapping Cord</u>: This component of the *MetaMoreFore* is a length of cord that is attached to a *Strapping-Specific Expansion Ring* or the *Strapping Tensioning Ring*.

Meatus: This is the opening in the glans through which urine flows.

Glans Cap: This is a component of the *MetaMoreFore* that covers the glans and is either used for "packing" of the *MetaMoreFore* or for use with the Bi-directional configuration. Early in R&D this was revered to as the Helmet (as it is reminiscent of the helmet that a Storm Trooper would wear).

<u>Push Piston</u>: This component of the **MetaMoreFore** is exclusivity used in a Bi-directional configuration. *Push Pistons* come in three lengths (10, 12 and 14 cm) and have notches along one side so as to latch the *Strapping Cord* of the *Strapping-Specific Expansion Ring*.